



AVIT
AARUPADAI VEEDU INSTITUTE OF TECHNOLOGY



VINAYAKA MISSION'S
RESEARCH FOUNDATION
(Deemed to be University under section 3 of the UGC Act 1956)



Accredited by NAAC



Approved by AICTE

Department of Biomedical Engineering

LAB MANUAL

Biomedical Instrumentation Lab

HOD(ECE/BME)

List of Experiments

1. Blood pressure measurement using sphygmomanometer
2. Design of instrumentation amplifier
3. Measurement PH using PH meter
4. Galvanic Skin resistance measurement
5. Recording of ECG using ECG simulator
6. Recording of EEG using EEG simulator
7. Recording of EMG using EMG simulator
8. Optical isolation Amplifier
9. Study of Phono Cardiogram (PCG)
10. Study of Types of electrodes

Ex. No 1 Blood Pressure Measurement

Aim: To Measure blood pressure using Sphygmomanometer and compare with digital blood pressure meter..

Apparatus Required:

1. Cuff
2. Inflator
3. Power supply
4. Stethoscope
5. Sphygmomanometer
6. Digital Blood pressure meter

THEORY:

Blood Pressure


Blood pressure is a measurement of the force applied to the walls of the arteries as the heart pumps blood through the body. The pressure is determined by the force and amount of blood pumped, and the size and flexibility of the arteries. Blood pressure is continually changing depending on activity, temperature, diet, emotional state, posture, physical state, and medication

The ventricles of heart have two states: systole (contraction) and diastole (relaxation). During diastole blood fills the ventricles and during systole the blood is pushed out of the heart into the arteries. The auricles contract anti-phase to the ventricles and chiefly serve to optimally fill the ventricles with blood. The corresponding pressure related to these states are referred to as systolic pressure and diastolic pressure .The range of systolic pressure can be from 90 mm of Hg to 145mm of Hg with the average being 120 mm of Hg. The diastolic pressure typically varies from 60mm of Hg to 90 mm of Hg and the average being 80 mmofHg.

1. SPHYGMOMANOMETER

Mercury Sphygmomanometer

This includes a mercury manometer, an upper arm cuff, a hand inflation bulb with a pressure control valve and requires the use of a stethoscope to listen to the Korotkoff sounds. Relies on the auscultatory technique.



2. Digital BP meter (Automated device)

This includes an electronic monitor with a pressure sensor, a digital display and an upper arm cuff. An electrically driven pump raises the pressure in the cuff. Devices may have a user-adjustable set inflation pressure or they will automatically inflate to the appropriate level, about 30 mmHg above the predicted systolic reading. On operation of the start button the device automatically inflates and deflates the cuff and displays the systolic and diastolic values. Pulse rate may so be displayed. Devices may also have a memory facility that stores the last measurement or up to 10 or more previous readings. It is battery powered and uses the oscillometric technique.



PROCEDURE

1. The upper arm is wrapped with the cuff belt connected to a mercury pressure gauge and air is pumped with a rubber ball to increase cuff pressure about 30 mmHg higher than the systolic blood pressure to block the artery and stop blood flow downstream.
2. Then, the cuff pressure is slowly lowered. The artery opens at the instant when the cuff pressure decreases below the systolic blood pressure and blood begins to flow on and off in synchrony with pulses causing the opening and closing of the artery. The sound emitted by the pulses is named Korotkoff's and continues until the cuff pressure decreases below the systolic blood pressure and the artery ceases the opening and closing.
3. The stethoscope placed closely to the artery downstream of the cuff is used to hear Korotkoff's sound; the blood pressures are measured. Cuff pressure when Korotkoff's sound begins to be heard is defined as the highest blood pressure (Systole) and that when the sound disappears is defined as the lowest pressure (Diastole).

Tabulation:

S.No	Patient Name	Sphygmomanometer		Semi Automated			Automated		
		Systolic (mmHg)	Diastolic (mmHg)	Systolic (mmHg)	Diastolic (mmHg)	Pulse (bpm)	Systolic (mmHg)	Diastolic (mmHg)	Pulse (bpm)
1	X								
2	Y								
3	Z								
4	A								

Mean Arterial Pressure (MAP)

$$MAP = (1/3 \times SBP) + (2/3 \times DBP) \quad \text{OR} \quad MAP = DBP + ((SBP - DBP)/3)$$

$$\text{Pulse pressure (PP)} = SBP - DBP$$

Hyper tension – Above normal BP; Hypo tension – Below normal BP

Result: Thus the blood pressure measurements are done using mercury sphygmomanometer and digital automated device.

Ex No. 2 Design of Instrumentation amplifier

AIM:

To construct, test, and perform instrumentation amplifier by using operational amplifier.


APPARATUS REQUIRED:

1. Op-amplifier Ic- 741
2. Resistor 1k ohm – 7
3. Dual RPS (0-30V)
4. Bread board
5. Connecting wires

THEORY:

An instrumentation amplifier is an integrated circuit (IC) that is used to amplify a signal. This type of amplifier is in the differential amplifier family because it amplifies the difference between two inputs. The importance of an instrumentation amplifier is that it can reduce unwanted noise that is picked up by the circuit. The ability to reject noise or unwanted signals common to all IC pins is called the common-mode rejection ratio (CMRR). Instrumentation amplifiers are very useful due to their high CMRR. Other characteristics, such as high open loop gain, low DC offset and low drift, make this IC very important in circuit design.

CIRCUIT DIAGRAM:



Gain calculation

$$V_o = (V_2 - V_1) A_v$$

$$\text{Where } A_v = \frac{R_4}{R_5} (1 + \frac{R_1 + R_3}{R_2})$$

PROCEDURE:

1. Connections are made as per circuit diagram
2. Check the connection correctly
3. Avoid loose connection and switch the power supply
4. Note down the output values for the given input and check the practical value with theoretical calculation

TABULATION

Sr. no	V1 volts	V2 volts	Theoretical Output Vo volts	Practical output Vo volts
1				
2				
3				
4				
5				

RESULT:

The instrumentation amplifier is designed and constructed and the output value was checked with theoretical calculation.

Ex NO. 3. Measurement of pH using pH meter

AIM

To measure the pH of a given sample

APPARATUS REQUIRED:


pH meter
pH Electrode
Duffer solution
Thermometer

THEORY:

Introduction: pH is a measure of the acidity of an aqueous solution. It is related to the concentration of hydrogen ion, H⁺. The pH scale can tell if a liquid is more acid or more base, just as the Fahrenheit or Celsius scale is used to measure temperature. The range of the pH scale is from 0 to 14 from very acidic to very basic. A pH of 7 is neutral. A pH less than 7 is acidic and greater than 7 is basic. Each whole pH value below 7 is ten times more acidic than the next higher value. The dissociation of water into hydrogen and hydroxide ions is the basis for the pH scale. The pH is related to the concentration of hydrogen ions by the formula, $\text{pH} = -\log[\text{H}^+]$. It is a logarithmic relationship. When the concentration of hydrogen ion changes by a factor of 10, the pH changes by a factor of 1. Many chemical processes are pH dependent and careful control of pH is important consideration. Solution pH is an important property that is measured by means of a pH meter.

A pH meter and its electrodes form a sensitive electrochemical device that will allow an accurate, reproducible, and reliable measurement of the pH of a solution (see Figure 1). A pH meter is essentially a voltmeter that measures the voltage of an electric current flowing through a solution between two electrodes. There is a direct relationship between the voltage and the pH of a solution. As a result, the meter on the instrument is calibrated directly in pH units. Two electrodes are required. One of them is called a glass electrode. This electrode is sensitive to the concentration of H⁺ ions in the solution. The other electrode is called the reference electrode and its operation is independent of the composition of the solution. The two electrodes are sometimes combined into a single entity called a combination electrode.

The pH Scale...



Formula:
 $\text{pH} = -\log(\text{H}^+)$



Figure 1. pH meter and a schematic diagram of a glass electrode.

PROCEDURE:

1. Prepare buffer solution of PH4 and Ph 9.2 by dissolving the respective buffer tablets
2. Measure the temperature of the buffer solution and adjust the temperature to penetrate.
3. Dip the electrode in the buffer solution.
4. Swich on the pH meter.
5. Push the Ph/Mv switch to ph positive and STBV/READ switch to read position and adjust CAL control to set on the read out and wait for 30 seconds
6. Set STBV/READ switch to STBV position to move the container with pH and buffer solution.
7. Wash the electrode with distilled water and clean with tissue paper.

PH MEASUREMENT OF GIVEN SAMPLE:

1. Measure the temperature of the given sample and set the temperature control and immersed electrode in the given sample.
2. Set the Ph/Mv switch to Ph position and STBV/READ switch to read position and wait for 30 seconds.
3. The pH value of the sample will be displayed in the read out.

Tabulation:

Sl.no	sample	pH
1	Sample 1	4
2	Sample 2	9.2
3	Sample 3	7

RESULT:

The PH value of the given solution was measured.

Ex No.4. Galvanic skin resistance measurements

AIM: To measure and display the galvanic skin resistance variations of skin.

APPARATUS REQUIRED:

1. GSR Simulator
2. GSR instrument
3. Electrode
4. DSO

PROCEDURE:

1. Connect the GSR meter to the main supply
2. Connect the resistance simulator to GSR meter
3. Set the resistance in simulator to 10 k ohms
4. Adjust the GSR meter to zero.
5. Increase the resistance value in simulator
6. Record the GSR meter and tabulated
7. Increase the resistance in simulator step wise and measure the deflection and tabulated.
8. Connect the electrodes in finger
9. Connect the output of electrodes to GSR meter
10. Measure the resistance change
11. Ask the patient for different emotions and record the change in GSR value.
12. Connect the DSO with GSR meter
13. Get the waveform of changes in GSR with respect to time.


THEORY:

Human skin offers resistance to electric current just like a resistor do in an electronic circuit. This transient change in the electrical conductivity of the skin followed by an arousal or oriented response is referred to as **Galvanic skin response or GSR**. The skin acts as a resistive layer and readily passes electric current. In normal Skin tone, the resistance varies from **25 Kilo Ohms to a few Mega Ohms**. This permits little current to pass through it. Sweating in the high emotional state or stress causes increased blood flow to the skin and reduces the skin resistance and increases its electrical conductivity.

Electrodermal activity (EDA), is the property of the human body that causes continuous variation in the electrical characteristics of the [skin](#). Historically, EDA has also been known as **skin conductance**, **galvanic skin response (GSR)**, **electrodermal response (EDR)**, **psychogalvanic reflex (PGR)**, **skin conductance response (SCR)**, **sympathetic skin response (SSR)** and **skin conductance level (SCL)**. The long history of research into the active and passive electrical properties of the skin by a variety of disciplines has resulted in an excess of names, now standardized to **electrodermal activity (EDA)**.^{[1][2][3]}

The traditional theory of EDA holds that skin resistance varies with the state of [sweat glands](#) in the skin. Sweating is controlled by the [sympathetic nervous system](#),^[4] and skin conductance is an indication of psychological or physiological [arousal](#). If the sympathetic branch of the [autonomic nervous system](#) is highly aroused, then sweat gland activity also increases, which in turn


increases skin conductance. In this way, skin conductance can be a measure of emotional and sympathetic responses




Procedure: Skin conductance (SC) is normally measured with 8mm diameter silver/silver chloride electrodes positioned on the medial phalanx of the index and middle fingers held in position by double sided sticky electrode collars. A non-saline jell is used. SC response (SCR) provides an indication of arousal.

Electrode placement for measurement of GSR is as shown in the fig. The output of electrode is given to the op amp input as shown in the circuit diagram.

The output of the amplifier if given to CRO to display the variations of skin resistance with respect to time. The sample waveform is as shown in the figure



Electrode placement diagram:



Tabular Colum

Basal Skin Resistance -

Sl.no	Simulator Skin resistance	Theoretical GSR	Measured GSR
1			
2			
3			
4			
5			

RESULT:

The galvanic skin resistance was recorded and analyzed the GSR for various emotions.

Ex. No 5: Generation and analysis of ECG wave using simulator

Aim:

To simulate ECG signal and analyze the signal parameters.

Apparatus Required:

ECG simulator, ECG Amplifier, Digital storage oscilloscope, connecting cables

Procedure

1. The ECG simulator is provided with 8 pattern selection
2. Connect the unit to mains.
3. Connect the simulator with ECG Amplifier
4. Connect the output of ECG amplifier to DSO.
5. Switch ON the Units.
6. Set the DSO setting to visualize the waveforms on DSO.
7. Note down the reading
8. Change the patterns to study different ECG waveforms

Theory:

The electrocardiogram is an instrument which records to electrical activity of heart. Electrical signals from the heart, characteristically placed the normal mechanical functions and monitoring of these signal has great clinical significance. Electrographs are used in authorization laboratories, coronary care units and for routine diagnostic applications in cardiology. The diagnostically useful frequency range is usually accepted as 0.05-150hz.



Einthoven Triangle:

It starts the vector sum of the projection of the frontal phase cardiac to vector to the 3 axis of the Einthoven triangle will be zero.

Lead 1: Left arm and right arm

Lead 2: Left leg and right leg

Lead3: Left and right arm



Bipolar Leads:

ECG is recorded by using two electrodes such that the final trace corresponds to the difference of electrical potentials existing between them, they are called standard leads and have been universally accepted. They are also called as einthoven lead.

Unipolar Lead:

If the electrode is placed close to the heart, higher potentials can be obtained, that normally available at limbs. ECG is recorded between a single exploratory electrode and the central terminal, which has

a potential corresponding to the centre of the body. The reference electrode is obtained by combination of several electrodes tied together at one point, it is of 2 types

- 1 Limb lead
- 2 Pericordial leads

LEAD 1

WAVE	AMP(V)	TIME(s)
P		
QRS		
T		

LEAD 2

WAVE	AMP(V)	TIME(s)
P		
QRS		
T		

Result:

The ECG was recorded and analyzed using simulator.

Ex. No. 6. Generation of EEG signals using EEG simulator.

Aim:

Generation of EEG signal using EEG simulator and measure the amplitude, frequency and to find the nature of the EEG.

APPARATUS REQUIRED:

1. EEG Simulator
2. EEG Amplifier
3. DSO
4. Connecting Cables

PROCEDURE

1. Connect the EEG Simulator to main
2. Connect the simulator output to eeg amplifier
3. Connect the EEG amplifier to DSO.
4. Switch ON the Units
5. Put the DSO on storage mode
6. Put the switch on DC mode
7. Time / Div Knob on the mS & S division
8. Voltage /Div Knob on the 5Mv
9. Now vary the amplitude as per requirement
- 10 If any noise on DSO, Check the ground

THEORY:

Electroencephalography (EEG) is the recording of electrical activity along the scalp produced by the firing of neurons within the brain. In conventional scalp EEG, the recording is obtained by placing electrodes on the scalp with a conductive gel or paste. Electrode locations and names are specified by the International 10–20 system for most clinical and research applications. Each electrode is connected to one input of a differential amplifier (one amplifier per pair of

electrodes); a common system reference electrode is connected to the other input of each differential amplifier. These amplifiers amplify the voltage between the active electrode and the reference.

A typical adult human EEG signal is about $10\mu\text{V}$ to $100\mu\text{V}$ in amplitude when measured from the scalp and is about $10\text{--}20\text{ mV}$ when measured from subdural electrodes.

EEG WAVE PATTERNS:





Fig. 13.5. Some examples of EEG waves.

DELTA WAVE:


Delta is the frequency range up to 4 Hz. It tends to be the highest in amplitude and the slowest waves. It is seen normally in adults in slow wave sleep. It is also seen normally in babies.



THETA WAVE:


Theta is the frequency range from 4 Hz to 7 Hz. Theta is seen normally in young children. It may

be seen in drowsiness or arousal in older children and adults; it can also be seen in meditation.



ALPHA WAVES:


Alpha is the frequency range from 8 Hz to 12 Hz. Hans Berger named the first rhythmic EEG activity he saw, the "alpha wave. It emerges with closing of the eyes and with relaxation, and attenuates with eye opening or mental exertion. The posterior basic rhythm is actually slower than 8 Hz in young children.



BETA WAVES:

Beta is the frequency range from 12 Hz to about 30 Hz Beta activity is closely linked to motor behaviour and is generally attenuated during active movements. It is the dominant rhythm in patients who are alert or anxious or who have their eyes open.

Since an EEG voltage signal represents a difference between the voltages at two electrodes, the display of the EEG for the reading encephalographer may be set up in one of several ways. The representation of the EEG channels is referred to as a montage.



Tabular Column

WAVES	AMPLITUDE(V)	TIME(S)	FREQUENCY(Hz)
ALPHA			
BETA			
THETA			
DELTA			

Result:

Thus the EEG waves are studied and the amplitude and time for each waveforms are noted for a Subject.

Ex No. 7. Recording of EMG using using EMG simulator

Aim

To record and analyze EMG signal using EMG simulator.

APPARATUS REQUIRED.

1. EMG simulator
2. Emg amplifier
3. DSO
4. Connecting Cables.

PROCEDURE

1. Connect the Biosignal amplifier to main supply.
2. Connect the EMG simulator to EMG amplifier
3. Switch on the Amplifier.
4. Connect the EMG amplifier out to DSO.
5. Switch ON the DSO.
6. Set Time/ Div knob ON THE mS division
7. Voltage division knob at 5Mv
8. Now vary the amplitude and frequency knob of EMG simulator
9. If any noise in signal check for proper ground.
10. Measure amplitude and frequency from the waveform.




Figure 2-2: Action potential (AP) of one motor unit

THEORY

The muscle cells are roughly cylindrical, with diameters between 10 and 100 μm but up to a few centimeters long. They may be arranged in parallel and bound by a connective tissue envelope into a homogeneous bundle. A myofiber is a multinucleated single muscle cell. It's basically water with some dissolved ions separated from the extra-cellular space that is mostly water with some dissolved ions. It generates a potential difference across its cell membrane by having different concentrations of ions. The fibers are excitable cells. Excitation signals are received at the synapse. Then a rapid depolarization occurs and is coupled with a contraction. It's a process during which electrochemical events occur. The action potentials are propagated along the sarcolemma, or cell membrane, toward the end of the fiber and downward from the surface into the transverse tubular system. The propagation of the action potential along a nerve or muscle fiber includes the flow of ions and gives rise to extra-cellularly recordable potential gradient. These potential gradients, moving in both time and space, constitute the electricity as recorded from active muscle fibers. Thus the small currents are generated prior to the generation of muscle force

Acquisition of EMG

As the brain's signal for contraction increases, it both recruits more motor units and increases the "firing frequency" of those units already recruited. All muscle cells within one motor unit become active at the same time. By varying the number of motor units that are active, the body can control the force of the muscle contraction. When individual motor contract, they repetitively emit a short burst of electrical activity known as the motor unit action potential (MUAP). It is detected by electrodes on the surface of the skin in proximity of the motor. The detection is illustrated in the following figure.

The function unit of a muscle is the motor. All the fibers which belong to one motor are activated at the same time. The motor unit action potential (MUAP) is the electrical response to the impulse from the axon. A MUAP looks like the following figure.




Figure 2-1: Detection of the motor unit action potential (MUAP)




Figure 2-2: Action potential (AP) of one motor unit

RESULT:

The EMG signal is recorded and analyzed the waveform.

ExNo.8. Opto-isolation amplifier

AIM


To measure pulse rate and to study the characteristics of opto isolation amplifier.

Apparatus required

1. Transducer : 1 with 1.5 meter, EP pin
2. Mains Cord : 1 Three pin with Socket
3. Test Probe : 2 2mm set.
4. Opto isolation Kit.
5. DSO
6. Connecting Cables

Installation Procedure

1. Connect the instrument to the mains.



2. Connect the Transducer to system

Insert the transducer plug into the transducer socket.



3. Use storage Oscilloscope.

1. Now switch ON DSO.

2. Put the DSO on storage mode,

3. Mode switch at DC position, Probe on 'X1' position.

4. Time / Div Knob on the 50mS Division. (Change the position as per signal)



5. Voltage / Div Knob on the 5V. (Change the position as per signal)

4. Now connect the DSO to the instrument.

Connect DSO probe to instrument at Pulse output 2mm red socket & Black socket Ground.

5. Putting the sensor

Put the finger on transducer & wrapped with Velcro around the finger. The transducer should not be too tight. Make arrangement of transducer as we get proper pulse. If there is any moments not get the proper readings. There are reading ± 2 digits up & down. Adjust threshold/Gain control if Pulse is not proper.



6. See signal Pattern on DSO

If Pulse not properly, adjust threshold control knob.

7. Display shows the Heart Rate. Red LED is indication of Heart Pulse. This same to Optical Isolation amplifier.


No movement, person must in rest position

8. Check the ground; if not proper make arrangement for that.

9. Now DSO Connect at isolation Output.

Green LED indication of Pulse output of Isolation amplifier. Red LED ON Green OFF & Pulse from Positive to Ground.

BLOCK DIAGRAM



MEASUREMENT OF PULSE RATE

Plethysmography is a non-invasive method for studies of the heart rate. This pulse wave will result in a change in the volume of arterial blood with each pulse beat. This change in blood volume can be detected in peripheral parts of the body such as the fingertip or ear lobe using a

technique called Photo plethysmography. The device that detects the signal is called a plethysmography.

Heart rate measurement indicates the soundness of the human cardiovascular system. The transducer consists of a Red LED that transmits an LDR through the fingertip of the subject, a part of which is reflected by the blood cells. The reflected signal is detected by a sensor. The changing blood volume with heartbeat results in a train of pulses at the output of the sensor, the magnitude is too small.

FORMULA

The Pulse to Pulse interval measures the period of heart beat is denoted by the following formula.

$$HR=60/T \text{ BPM}$$

T=Pulse to pulse Interval in Second. (P-P)

Where HR is the heart rate measured in beat-per-minute (BPM), interval measured in millisecond (ms). For example, if P-P is 800 ms, the heart rate is 75 BPM.

TABULAR COLUMN

Sl.No	Subject	Pulse time ms	Calculated HR	HR Display
1				
2				
3				
4				
5				

THEORY


Opto-isolators or Opto-couplers, are made up of a light emitting device, and a light sensitive device, all wrapped up in one package, but with no electrical connection between the two, just a beam of light. The light emitter is nearly always an LED. The light sensitive device may be a photodiode, phototransistor, or more esoteric devices such as thyristors, TRIACs etc.

The opto-coupler is a component that contains the two elements required for an opto-isolator:

Light emitter: The light emitter is on the input side and takes the incoming signal and converts it into a light signal. Typically the light emitter is a light emitting diode.

Light detector: The light detector within the opto-coupler or opto-isolator detects the light from the emitter and converts it back into an electrical signal. The light detector can be any one of a number of different types of device from a photodiode to a phototransistor, photodarlington, etc.

Opto-coupler symbol:




The symbol shows the LED, which is normally used as the light emitter. The opto-coupler symbol also shows the receiver, often a phototransistor or photodarlington, although other devices including light sensitive diacs, etc may also be used. The relevant device type is shown within the optocoupler circuit symbol.

Application of MCT2E:

It is a combination of 1 LED and a transistor. Pin 6 of transistor is not generally used and when light falls on the base-emitter junction then it switches and pin5 goes to zero.

MCT2E Opto-Coupler - Circuit



When logic zero is given as input then the light doesn't fall on transistor so it doesn't conduct which gives logic one as output.

When logic 1 is given as input then light falls on transistor so that it conducts, that makes transistor switched ON and it forms short circuit this makes the output is logic zero as collector of transistor is connected to ground.

RESULT

The Characteristics of opto-isolation amplifier was studied and Pulse rate was measured.

Ex. No 9 Study of Phono Cardiogram for measurement of Heart sounds

AIM

The basic aim of phonocardiograph is to pick up the different heart sounds, filter out the heart sounds and to display or record them.


THEORY

A Phonocardiogram or PCG is a plot of high fidelity recording of the sounds and murmurs made by the heart with the help of the machine called phonocardiograph, or "Recording of the sounds made by the heart during a cardiac cycle". The sounds are thought to result from vibrations created by closure of the heart valves. There are at least two: the first when the atrioventricular valves close at the beginning of systole and the second when the aortic valve closes at the end of systole. It allows the detection of sub audible sounds and murmurs, and makes a permanent record of these events. In contrast, the ordinary stethoscope cannot detect such sounds or murmurs, and provides no record of their occurrence. The ability to quantitate the sounds made by the heart provides information not readily available from more sophisticated tests, and provides vital information about the effects of certain cardiac drugs upon the heart. It is also an effective method for tracking the progress of the patient's disease.


Heart sounds are classified into four groups on the basis of their mechanism of origin, they are

1. Valve closure sound
2. Ventricular filling sound
3. Valve opening sounds and
4. Extra cardiac sounds

HEART SOUNDS



BLOCK DIAGRAM



These sounds occur at the beginning of systole (first heart sound) and the beginning of diastole (second heart sound). The first heart sound is due to the closure of mitral and tricuspid valves associated with myocardial contraction. And the second heart sounds is due to the closure of the aortic and pulmonary valves. The first heart sounds are low frequency vibrations occur approximately 0.05s after the onset of the QRS complex of the ECG, the first heart sounds last for (0.1 to 0.12s) and the frequency ranges 30-50Hz. The second heart sound is due to the vibrations set up by the closure of semilunar valves. These sounds start approximately (0.03 to 0.05)s after the end of T wave of the ECG, this lasts for (0.08 to 0.14)s and have a frequency up to 250Hz.

Ventricular filling sounds

These sounds occur either at the period of rapid filling of the ventricles (third heart sound) or during the terminal phase of ventricular filling. These sounds are inaudible. Third heart sound starts at (0.12 to 0.18) s after the onset of the second heart sound. it last approximately (0.04 to 0.08) s. The frequency is about 10 to 100 Hz.

Valve opening sounds

These sounds occur at the time of opening of the atria ventricular valves and semi lunar valves. The fourth heart sound starts approximately (0.12 to 0.18) s after the onset of the P wave. The sound last for (0.03 to 0.06)s. And the frequency is 10 to 50 Hz.

Extra cardiac sounds

These sounds occur in late systole or early diastole and are believed to be caused by thickened pericardium which limits ventricular distensibility. Murmurs are sounds related to non-laminar flow of blood in the heart and the great vessels. They are distinguished from the basic heart sounds such that they have noisy character having long duration and with high frequency components up to 1000 Hz.

RESULT:

Thus, In this study, by using a Phonocardiograph, different heart sounds have been identified and analyzed.

Ex No. 10. Study of Different types of Electrodes,

Aim: To study the different types of electrodes.

Theory:

Biopotential electrodes is a transducer that convert the body ionic current in the body into the traditional electronic current flowing in the electrode. Biopotential electrode should be able to conduct small current across the interface between the body and the electronic measuring circuit. We will get to know Electrical characteristics of biopotential electrodes, Different type of biopotential electrodes, Electrodes used for ECG, EEG, EMG, and intracellular electrodes.

1. Microelectrodes: Electrodes use to measure bioelectric potential near or within a single cell .
2. Skin surface electrodes: Electrodes used to measure ECG , EMG , EEG potential from the surface of the skin.
3. Needle electrodes: Electrodes are used to penetrate the skin to record EEG potential from a local region of the brain or EMG potential from gray.

MICRO ELECTRODES:

There are electrodes with tip sufficiently small to penetrate a single cell in order to obtain reading from the cell . They have tip diameters ranging from approx 0.05 to 10 microns. They can be formed from solid metal needles.

Types of Microelectrodes:

Metal microelectrodes:

They are formed by electrolytic etching the tip of a fine tungsten or stainless steel wire with a desired size then the wire is coated at most to the tip with an insulating material.

The structure of a metal microelectrode for intracellular recordings.




Figure 5.18 Structures of two supported metal microelectrodes (a) Metal-filled glass micropipet. (b) Glass micropipet or probe, coated with metal film.


Fig. 1: Metal Microelectrode

Micropipette microelectrode:

This type is a glass with tip drawn out to desired size it is fitted with an electrolyte compatible with cellular fluids . It has dual interface.

1. Metal wire in contact with electrolyte solution inside micropipette.
2. Electrode inside the pipette and fluids inside or immediately outside the cell.

Microelectrodes



A glass micropipet electrode filled with an electrolytic solution (a) Section of fine-bore glass capillary. (b) Capillary narrowed through heating and stretching. (c) Final structure of glass-pipet microelectrode.

BODY SURFACE ELECTRODES:

Electrodes used to obtain bioelectric potential from the surface of the body are called as body surface electrodes.

Types of body surface electrodes:

I. Metal plate electrodes:

It consist of a flat metal plate that has bent into a cylindrical segment. A terminal is placed on its outside surface near 1 end. This terminal is attached to the lead wire to the EEG.

II. Metal disk electrode:


This electrode which has a lead wire welded to the back surface made up of several different material. It is coated with electrolyte gel and then passed against patients chest wall.

III. Disposable foam pad electrode:

It consist of relatively large disk of plastic foam material with a silver disk on one side attached to the silver plated snap similar to that of lead wire with female portion is then snapped on electrode.

Body-Surface Recording Electrode

Metal-Plate Electrodes




Body-surface biopotential electrodes (a) Metal-plate electrode used for application to limbs. (b) Metal-disk electrode applied with surgical tape. (c) Disposable foam-pad electrodes, often used with electrocardiograph monitoring apparatus.

IV. Suction cup electrode:

This electrodes are frequently used in ECG at pre-cardiac levels because they can be placed at particular location and used to take a recording that consist of a hollow metallic cylindrical electrode that makes contact with skin at its base. In appropriate technical form the lead wire is attached to metal cylinder and rubber suction bulb it fits over its oyer base. Electrolyte gel is placed over the contacting surface of electrodes.

Body-Surface Recording Electrode

Suction Electrodes




A metallic suction electrode is often used as a precordial electrode on clinical electrocardiographs. **No need for strap or adhesive** and can be **used frequently**. **Higher source impedance** since the contact area is small

FLOATING ELECTRODE:

The principle feature of this electrode is that the actual electrode element is released in cavity so that it does not come in contact with skin itself. Instead the element is surrounded by electrolyte gel in cavity.

Floating Electrodes



The recess in this electrode is formed from an open foam disk, saturated with electrolyte gel and placed over the metal electrode. **Minimize motion artifact**

INSULATED NEEDLE ELECTRODE:

The basic needle electrode consists of a solid needle usually made up of stainless steel with a sharp point insulated with a coating such as insulating varnish and its tip is exposed.

Types of needle electrodes:

I. Co-axial needle electrode:

It consists of a small gauge hypodermic needle that is modified by running an insulated fine wire down the centre of its lumen and filling the remainder of its lumen with an electrode material.

II. Bipolar co-axial electrode:

The two wires are placed within the lumen of a needle and can be connected differentially so as to be sensitive to electrical activity only in the immediate vicinity of the electrode tip.




Fig. 6: Needle Electrodes

Questions:

- What are the different types of bio-potential electrodes?
- Why are the types of Microelectrodes?
- What are the applications of skin surface electrodes?

Result:

The types electrodes and their uses are studied